Fast, safe, non-contact measurement of human body temperature

- Display in Temperature in °C & °F
- Forehead Non-Contact
- Range 32°C -42°C
- CE, RoHS

Simply point the non-contact IR (infra-red) thermometer at the subject's forehead from a distance of 2" to 6", pull the trigger and, within 1 second, obtain immediate, highly accurate readings of human body temperature. Automatically adjusts for the difference between the surface temperature of the forehead and internal body temperature of the subject, for highly accurate reading within +/-0.2°C.



No. Description

DT-8826 Clinical IR Thermometer

	Range	Resolution	Accuracy
Temperature °F	89.6 to 106ºF	0.1°F	±0.4ºF
Temperature °C	32 to 42.0ºC	0.1°C	±0.2°C

